

Fully Edentulous Maxilla or Mandible

Step 1: On a dry, well fitting denture, place 6 x-ray markers.

- -Place the markers on varying axial planes.
- -Example: If you place a marker on the buccal area of #2-3, then place your next marker on the lingual area of #4-5 area, the next marker on the facial of #8-9 area and so on.

Do not place the markers on the intaglio surface of the denture







Step 2- Pre-Scan Patient Instructions:

-Use a bite registration for separation, secure fit, and to eliminate rocking
-Confirm fit by asking the patient, "Is the denture fitting correctly?"
-Breathe through your nose, close your eyes, and place the tongue at the back of the mouth

Step 3- First Scan: Capture CBCT of ONLY the maxilla denture with the x-ray markers

-Position appliance on a foam block -Orientate teeth occlusal DOWN CBCT General Settings: 6cm-8cm, 0.3 – 0.4 Voxel

Step 4. Second Scan: Capture CBCT of the patient WEARING maxilla denture with x-ray markers.

Vital Anatomy to be captured:

-Superior: Mid Sinus / below orbit

-Inferior: Capture 1-3 mm of mandible dentition

CBCT General Settings: Maxl 6cm-8cm, 0.3-0.4 Voxel

Step 5. When scanning a patient following Dual Scan Protocol:

-Review the scan carefully to ensure there is no "Black Air Space" between the denture and the soft tissue of the patient



Step 6. Export and Upload CBCT Scans to Implant Concierge

-Export the CBCT scans individually, in a multi-file DICOM format to the desktop

- -Right click on each folder, go to "Send To", and select "Compressed (Zipped) Folder"
- -Upload both zipped/compressed folders to Implant Concierge



Placement of the Gutta Percha Markers

- Ensure the denture surface is dry before placing the sticker markers
- Use a well-fitting, hard lined denture with NO soft liner; your surgical guide will only fit as well as the denture does without soft liner!
- Place at least 6 markers (no need to exceed 8 markers) on varying planes split between the buccal and palatal/lingual, making sure they do not overlap with one another. For example: if you place a marker on the buccal area of #2 then place your next marker on the palatal side of #4/5 area, then place the next marker on the facial of #8/9 area and so on. Do not place the markers on the intaglio surface of the denture.
- Place the metal BB just above the gingival margin of the denture teeth. Do not place markers on the palatal area.
- Remember you will be doing two scans: One scan of the marked denture by itself on a non-radiopaque surface like packing foam, and another of the patient wearing the marked denture. Do not remove the markers in between scans.





Partially Edentulous Maxilla & Mandible

Step 1. The treating doctor performs initial physical examination of the patient and takes a digital impression of both arches.

Step 2. Patient Positioning:

- Breathe through your nose
- Bite gently on cotton rolls
- Superior: Sinus (below orbit)
- Inferior: Chin and entire mandible



Step 3. CBCT General Settings: Full 13cm-20cm, 0.3-0.4 Voxel

Step 4. Export and Upload CBCT Scan

-Export the CBCT scan in a multi-file DICOM format to the desktop

-Right click on the folder, go to "Send To", and select "Compressed (Zipped) Folder"

-Upload the zipped/compressed folder to Implant Concierge

Step 5. Export & Upload STL File of the Digital Impression