


# LAB SCREW Data Sheet

ZERAMEX® **P6**



## SN Platform




### **P35450**

**Matches:**  
P6 Abutment SN straight, 1mm (P15501)  
P6 Abutment SN straight, 2mm (P15502)

**Distinguishing feature to other SN screws:**  
Length: 7.4mm

**Tightening torque:** 5 Ncm



### **P35451**


**Matches:**  
P6 Abutment SN angular 15° (P15515)

**Distinguishing feature to other SN screws:**  
Length : 6.3mm

**Tightening torque:** 5 Ncm



## RN Platform




### **P36450**

**Matches:**  
P6 Abutment RN straight, 1mm (P16501)  
P6 Abutment RN straight, 2mm (P16502)

**Distinguishing feature to other RN screws:**  
Length : 7.5mm

**Tightening torque:** 5 Ncm



### **P36451**

**Matches:**  
P6 Abutment RN angular 15° (P16515)

**Distinguishing feature to other RN screws:**  
Length : 6.7mm

**Tightening torque:** 5 Ncm