




ZERAMEX[®] (P6)

Product range

PRODUCT RANGE

ZERAMEX®P6 Implants

Ø 3.3 mm	Art. No.	Name	Dimension	Material
	P15508	ZERAMEX®P6 Implant Ø 3.3 x 8 mm SN (incl. healing cap)	Length: 8mm	ZrO ₂ -ATZ-HIP
	P15510	ZERAMEX®P6 Implant Ø 3.3 x 10 mm SN (incl. healing cap)	Length: 10mm	
	P15512	ZERAMEX®P6 Implant Ø 3.3 x 12 mm SN (incl. healing cap)	Length: 12mm	
Ø 4.1 mm	Art. No.	Name	Dimension	Material
	P16508	ZERAMEX®P6 Implant Ø 4.1 x 8 mm RN (incl. healing cap)	Length: 8mm	ZrO ₂ -ATZ-HIP
	P16510	ZERAMEX®P6 Implant Ø 4.1 x 10 mm RN (incl. healing cap)	Length: 10mm	
	P16512	ZERAMEX®P6 Implant Ø 4.1 x 12 mm RN (incl. healing cap)	Length: 12mm	
Ø 4.8 mm	Art. No.	Name	Dimension	Material
	P17508	ZERAMEX®P6 Implant Ø 4.8 x 8 mm RN (incl. healing cap)	Length: 8mm	ZrO ₂ -ATZ-HIP
	P17510	ZERAMEX®P6 Implant Ø 4.8 x 10 mm RN (incl. healing cap)	Length: 10mm	
	P17512	ZERAMEX®P6 Implant Ø 4.8 x 12 mm RN (incl. healing cap)	Length: 12mm	



ZERAMEX® P6 Abutments

SN Platform	Art. No.	Name	Dimension	Material
	P15501	ZERAMEX®P6 Abutment SN Straight, 1 mm (incl. screw)	AH: 6.5mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P15502	ZERAMEX®P6 Abutment SN Straight, 2 mm (incl. screw)	AH: 6.5mm, CH: 2mm	
	P15515	ZERAMEX®P6 Abutment SN Angular 15°, 1 mm (incl. screw)	AH: 5.5mm	ZrO ₂ -ATZ-HIP VICARBO®


CAD/CAM

	P15530	ZERAMEX®P6 ZERABASE SN, for crown (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P15531	ZERAMEX®P6 ZERABASE SN, for bridge (incl. screw)	AH: 7mm, CH: 1mm	

RN Platform	Art. No.	Name	Dimension	Material
	P16501	ZERAMEX®P6 Abutment RN Straight, 1 mm (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P16502	ZERAMEX®P6 Abutment RN Straight, 2 mm (incl. screw)	AH: 8mm, CH: 2mm	
	P16515	ZERAMEX®P6 Abutment RN Angular 15°, 1 mm (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®

CAD/CAM

	P16530	ZERAMEX®P6 ZERABASE RN, for crown (incl. screw)	AH: 7mm, CH: 1mm	ZrO ₂ -ATZ-HIP VICARBO®
	P16531	ZERAMEX®P6 ZERABASE RN, for bridge (incl. screw)	AH: 7mm, CH: 1mm	



VICARBO® screws	Art. No.	Name	Dimension	Material
	P15450	ZERAMEX®P6 VICARBO SN Straight	Length: 7.3mm	VICARBO®
	P15451	ZERAMEX®P6 VICARBO SN Angular	Length: 6.25mm	
	P16450	ZERAMEX®P6 VICARBO RN Straight	Length: 7.5mm	
	P16451	ZERAMEX®P6 VICARBO RN Angular	Length: 6.7mm	

ZERAMEX® ZERALOCK™ Abutments

Straight & angular abutments

SN Platform	Art. No.	Name	Dimension	Material
	P15305	ZERAMEX® ZERALOCK™ Abutment SN straight	AH: 6mm	ZrO ₂ -ATZ-HIP
	T15315	ZERAMEX® ZERALOCK™ Abutment SN angular 15°	AH: 5.4mm	

RN Platform	Art. No.	Name	Dimension	Material
	P16305	ZERAMEX® ZERALOCK™ Abutment RN straight	AH: 6mm	ZrO ₂ -ATZ-HIP
	T16315	ZERAMEX® ZERALOCK™ Abutment RN angular 15°	AH: 5.5mm	



LOCATOR®

SN Platform	Art. No.	Name	Dimension	Material
	P15202	ZERAMEX® ZERALOCK™ Locator® 2 mm SN	Height: 2mm	ZrO ₂ -ATZ-HIP
	P15203	ZERAMEX® ZERALOCK™ Locator® 3 mm SN	Height: 3mm	
	P15205	ZERAMEX® ZERALOCK™ Locator® 5 mm SN	Height: 5mm	

RN Platform	Art. No.	Name	Dimension	Material
	P16202	ZERAMEX® ZERALOCK™ Locator® 2 mm RN	Height: 2mm	ZrO ₂ -ATZ-HIP
	P16203	ZERAMEX® ZERALOCK™ Locator® 3 mm RN	Height: 3mm	
	P16205	ZERAMEX® ZERALOCK™ Locator® 5 mm RN	Height: 5mm	



ZERAMEX® Prosthetics - SN Platform

Soft Tissue Management	Art. No.	Name	Dimension	Material
	P35500	ZERAMEX®P6 Healing Cap SN	Height: 0.5mm	PEEK
	P35503	ZERAMEX®P6 Gingiva Former SN, 3 mm	Height: 3mm	
	P35504	ZERAMEX®P6 Gingiva Former SN, 4 mm	Height: 4mm	
	P35530	ZERAMEX®P6 Provisional SN (incl. screw)	AH: 7mm CH: 1mm	PEEK, PEEK-CW30

Impression Taking	Art. No.	Name	Dimension	Material
	P35510	ZERAMEX®P6 Transfer Open Tray SN	Sleeve height: 10mm Height incl. screw: 19mm	PEEK-CW30, Aluminum
	P35512	ZERAMEX®P6 Transfer Closed Tray SN	Sleeve height: 7mm Height incl. screw: 14mm	PEEK-CW30, Aluminum
	P35513	ZERAMEX®P6 Transfer Closed Tray, long SN	Sleeve height: 11mm Height incl. screw: 19mm	PEEK-CW30, Aluminum
	P35514	ZERAMEX®P6 Scanbody SN	Height: 10mm	PEEK, PEEK-CW30

PRODUCT RANGE

ZERAMEX® Prosthetics - Platform

Auxiliaries	Art. No.	Name	Dimension	Material
	P35520	ZERAMEX®P6 Implant Replica SN	Height: 10mm	Aluminum
	P35550	ZERAMEX®P6 Provisional Screw SN	Length: 7.3mm	PEEK-CW30
	P35553	ZERAMEX®P6 Lab Screw SN	Length: 6.2mm	PEEK

ZERAMEX® Prosthetics - Platform

Soft Tissue Management	Art. No.	Name	Dimension	Material
	P36500	ZERAMEX®P6 Healing Cap RN	Height: 0.6mm	PEEK
	P36502	ZERAMEX®P6 Closure Screw RN	Height: 0.5mm	
	P36503	ZERAMEX®P6 Gingiva Former RN, 3 mm	Height: 3mm	
	P36504	ZERAMEX®P6 Gingiva Former RN, 4 mm	Height: 4mm	PEEK, PEEK-CW30
	P36530	ZERAMEX®P6 Provisional RN (incl. screw)	AH: 7mm CH: 1mm	

Impression Taking	Art. No.	Name	Dimension	Material
	P36510	ZERAMEX®P6 Transfer Open Tray RN	Sleeve height: 11mm Height incl. screw: 20mm	PEEK-CW30, Aluminum
	P36512	ZERAMEX®P6 Transfer Closed Tray RN	Sleeve height: 7mm Height incl. screw: 14mm	PEEK-CW30, Aluminum
	P36513	ZERAMEX®P6 Transfer Closed Tray, long RN	Sleeve height: 11mm Height incl. screw: 18mm	PEEK-CW30, Aluminum
	P36514	ZERAMEX®P6 Scanbody RN	Height: 10mm	PEEK, PEEK-CW30


Auxiliaries	Art. No.	Name	Dimension	Material
	P36520	ZERAMEX®P6 Implant Replica RN	Height: 10mm	Aluminum
	P36550	ZERAMEX®P6 Provisional Screw RN	Length: 6.1mm	PEEK-CW30
	P36553	ZERAMEX®P6 Lab Screw RN	Length: 5.2mm	PEEK
	P18550	ZERAMEX®P6 Abutment Planning Kit		PEEK

ZERAMEX® Tools

	Art. No.	Name	Material
	P35601	ZERAMEX®P Rose Drill ø 2.2 mm	Stainless steel
	P35602	ZERAMEX®P Pilot Drill ø 2.2 mm	
	P35633	ZERAMEX®P Drill ø 2.8 mm	
	P36633	ZERAMEX®P Drill ø 3.5 mm	
	P37633	ZERAMEX®P Drill ø 4.2 mm	
	P35622	ZERAMEX®P Profile Drill SN ø 2.8 mm	
	P36622	ZERAMEX®P Profile Drill ø 3.5 mm	
	P35650	ZERAMEX®P Depth Gauge ø 2.0 mm	
	P36650	ZERAMEX®P Depth Gauge ø 3.0 mm	
	P35620	ZERAMEX®P Screw Tap ø 3.3 mm	
	P36620	ZERAMEX®P Screw Tap ø 4.1 mm	Stainless steel
	P37620	ZERAMEX®P Screw Tap ø 4.8 mm	
	P38619	ZERAMEX® Prosthetic Wrench, 19 mm	
	P38623	ZERAMEX® Prosthetic Wrench, 23 mm	
	T38620	ZERAMEX® Pickup 20 mm	
	T38622	ZERAMEX® Rescue Pickup	
	T38625	ZERAMEX® Pickup 25 mm	
	T38626	ZERAMEX® Alignment Tool	
	C7650	Drill extension	

PRODUCT RANGE

ZERAMEX® Tools

	Art. No.	Name	Material
	P35453	ZERAFIX™ Placeholder SN 3.0mm	PEEK
	P36453	ZERAFIX™ Placeholder RN 3.3mm	
	P38645S	ZERALOCK™ Abutment Holder Set SN & RN	PEEK
	P38644S	ZERAFIX™ Abutment Holder Set SN & RN	
	P48850	ZERAMEX®P Surgery Tray	Stainless steel
	P48854	ZERAMEX®P Surgery Tray, fully equipped with ratchet	
	P48870	ZERAMEX®P Surgery Tray "mini"	
	P48932	ZERAMEX®P Ratchet Adapter Unit Short	
	P48935	ZERAMEX®P Surgical Ratchet	
	P48860	ZERAMEX®P Prosthetics Tray	
	P48865	ZERAMEX®P Prosthetics Kit, incl. ratchet	

ZERAMEX® Locator® Accessories

	Art.No.	Name	Material
	XT38253	Docklocs® laboratory set, up to 40° divergence compensation: 2 zirconium housings (Ø 5.8 mm, height 2.5 mm) with black process insert (height 1.9 mm), 2 pieces block-out ring, 2 pieces replacement males, green, 2 pieces replacement males, orange, 2 pieces replacement males, red	Santoprene® TPE/Polyamid/Housing ZiO2 HD-PE Purell
	XT38251	Docklocs® laboratory set, up to 40° divergence compensation: 2 metal housings (Ø 5.5 mm, height 2.5 mm) with black process insert (height 1.9 mm), 2 pieces block-out ring, 2 pieces replacement males, green, 2 pieces replacement males, orange, 2 pieces replacement males, red	Santoprene® TPE/Polyamid/Housing Titan G5 HD-PE Purell
	XT38205	Docklocs® Replacement male, red, extra-light retention, 10°–20°, 8 pcs.	Polyamid
	XT38206	Docklocs® Replacement male, orange, light retention, 10°–20°, 8 pcs.	Polyamid
	XT38207	Docklocs® Replacement male, green, strong retention, 10°–20°, 8 pcs.	Polyamid
	XT38209	Docklocs® Blockout Ring, 20 pcs.	Santoprene® TPE
	XT38230	Zirconium housing with processing insert, 2 pcs.	Housing ZiO2 HD-PE Purell
	XT38210	Metal housing titanium with processing insert, 4 pcs.	Housing Titan G5 HD-PE Purell
	XT38214	Docklocs® lab analog straight (Ø 4 mm), 4 pcs.	Titan Grade 5
	XT38215	Docklocs® impression coping with black processing insert, 4 pcs.	Housing Titan G5 HD-PE Purell

CE	0050	⚠	Note (observe accompanying documents)
REF	Order number	📅	Expiry date:
LOT	Batch ID	📖	Follow the instructions for use
🔄	Do not reuse	🏭	Manufacturing date:
STERILE	Steam sterilized	🏢	Manufacturer
🧼	Unsterile	🚫	Do not use if packaging is damaged
☀️	Keep away from sunlight	FDA	The products of the portfolio are FDA approved
🌧️	Keep dry	Rx only	Prescription use only
🇺🇸	Implant small ø3,3mm SN (Platform ø4,0mm)	🇧🇪	Implant small ø4,1mm RN (Platform ø4,8mm)
🇪🇺	Implant small ø4,8mm RN (Platform ø4,8mm)		